



**SAMPILIO X8™**

Dissolution Autosampler

EFFICIENT · ACCURATE · RELIABLE · MULTIPURPOSE

## SAMPILIO X8

Dissolution Autosampler

EFFICIENT · ACCURATE · RELIABLE · MULTIPURPOSE



### INTRODUCTION

Normally comprises a sample collector, a peristaltic or syringe pump to provide the motive force to transport the samples from the dissolution tester to the collector, a PC and interface box to control the system during operation.

The principle of operation is simple - medium from each of the dissolution vessels is circulated via an 6- or 8-line peristaltic pump through 6 or 8 switching valves prior to being returned to the dissolution vessel.

At user-defined intervals the valves operate, diverting a preset volume of sample into the sample collection lines, whereupon the samples are dispensed into either test tubes or open HPLC vials (or injected directly into sealed septum vials by means of an electrically operated vial piercing head provided for that purpose).

The pump is then reversed to clear the sampling lines prior to the next sampling interval, whereupon the operation is repeated. The whole operation is controlled and moni-

tored by a PC. The exact status of the test at any given time can be determined from the software.

In the case of test tubes, the samples must be handled manually, for example by presenting them to the „sipper“ accessory of a suitable spectrophotometer.


HPLC vials containing samples can be removed at any time and placed directly into an HPLC Autosampler.

This version is particularly useful where analytical techniques other than UV/Vis or HPLC are employed or where the samples require a degree of manipulation, for example, to be diluted or mixed with a reagent prior to analysis.

### PRINCIPLE OF OPERATION

The „Sampilio X8“ is a new generation of „Off-Line“ Dissolution Sampling System specifically designed with a series of 6, 7 or 8 dedicated bidirectional small volume diaphragm pumps (one per line/vessel) to facilitate the flush-sample-purge functions.





As well as being extremely accurate (Volumetric Precision <math><0.25\text{ mL}</math>, typically

- First In/First Out (FIFO) principle
- Low dead volume
- Eliminates need for media replacement
- Low cross contamination
- Short sampling interval times (2 min)
- long life for years with no abrasion
- wide range of applications, incl. organic solvents usage

The First In/First Out (FIFO) principle employed in the system is the same as that found in manual testing.

The low dead volumes employed in the system ensure that flush, sample and purge times are kept to a minimum whilst flush media recycling makes filter changer and media replacement obsolete and dissolution calculations simple. Cross contamination is <math><1\%</math> at 2 minute sampling intervals. The short

interval time is particularly important when testing quick release formulations in so much that it allows sampling at intervals hitherto unachievable by more conventional methods. The user interface is simple, functional and easy to use.

The unit is supplied as standard with two collection racks, one to accept 2 mL HPLC vials and the other 10 mL test tubes.

Each rack accommodates 10 rows of 8 lines and an additional row with test tubes for waste. With additional racks (optional) it is possible to increase up to 100 samples each line.

In order to eliminate any cross contamination, the standard sampling procedure is always flush-sample-purge, so that all not sampled media is completely recycled.

The unique delivery system with individual bidirectional pumps offers a wide range of applications like staggered sampling or transfer to UV, HPLC (option).

The option of cooling the collected sample makes the Sampilio X8 the most usable Autosampler available.

## SPECIFICATIONS

<b>DIMENSIONS</b>	B 30cm x H 55cm x D 58cm
<b>WEIGHT</b>	23kg tare weight ( 28 kg packing included)
<b>VOLTAGE SUPPLY AND POWER RATING</b>	90 - 250 V 50/60 Hz, 250 W
<b>ENVIRONMENTAL CONDITIONS</b>	15 - 40°C, less 95% humidity
<b>INTERFACES</b>	USB, LAN, RS-232, RS-485, Contact Closure (WAGO)
<b>INTERNAL METHOD STORAGE</b>	99 methods
<b>TEST TUBE AND VIAL CAPACITY</b>	HPLC vials 2mL (max. OD 12mm, max H 35mm), test tubes 10mL (max OD 14mm, max H 100mm)
<b>RACKS</b>	Test Tube Racks and Vial Racks with 10 rows and 8 lines, 1 additional line for Waste with Test Tubes, material: POM
<b>RACK CHANGE</b>	possible during running method
<b>PUMP FLOW RATE</b>	1 - 25 mL/min (limited by used ID of tubings: 1.0 or 1.2 mm version available)
<b>VOLUMETRIC PRECISION</b>	<0.25 mL, typical 0.1 mL
<b>CROSS CONTAMINATION</b>	< 1 % at 2 min sampling interval
<b>SAMPLE INTERVAL TIME</b>	2min - 999h (limited on selected flow rates)
<b>REMOTE CONTROL SOFTWARE</b>	no installation required, control via Web browser
<b>TUBINGS VOLUME</b>	without pumps, internal tubings, sampling needles and filters, the volumes are: Low Flow Tubing Kit with 1.0 mm ID: 1.0mL High Flow Tubing Kit with 1.2mm ID: 2.26mL
<b>INTERNAL VOLUME</b>	pumps, needles and internal tubings volumes are: 0.5mL
<b>TUBING TYPE</b>	PTFE, 1/16" OD with bend protection
<b>TEMPERATURE OPTION</b>	Peltier technique with built in power supply, 5 - 37°C (limited on environment conditions)
<b>TRANSFER OPTION</b>	Option to transfer samples to UV or HPLC with an internal Transfer Module (triggering from UV or HPLC possible)